

BornAgain - Bug #2076

Beam distribution editors are not updating on parameter change

13 Jun 2018 16:38 - pospelov

Status:	Resolved	Start date:	13 Jun 2018
Priority:	Normal	Due date:	
Assignee:	pospelov	% Done:	0%
Category:		Estimated time:	0.00 hour
Target version:	Sprint 38		

Description

- Go to instrument view, call "Extended distribution viewer" for beam wavelength.
- Switch from monochromatic to Gaussian distribution
- Step inside "Number of samples" integer spin box and try to change number of samples
- Distribution still show number of samples 5

To make it show new value one have to switch to another widget

Before it was not like that. Changing the value in combo box was immediately propagated to custom plot.

History

#1 - 13 Jun 2018 16:38 - pospelov

- Subject changed from *Beam distribution editors do not update on parameter change* to *Beam distribution editors are not updating on parameter change*

#2 - 19 Jun 2018 09:17 - pospelov

- Assignee set to pospelov

#3 - 19 Jun 2018 16:57 - pospelov

- Status changed from *Sprint* to *Resolved*